

## Product datasheet

# Sheep Anti-Cow IgG heavy chain (HRP) ab112618

### Overview

Product name	Sheep Anti-Cow IgG heavy chain (HRP)
Host species	Sheep
Target species	Cow
Tested applications	<b>Suitable for:</b> WB, IHC-P, ICC/IF, ELISA
Conjugation	HRP

### Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 6.8 Constituent: 0.05% CMIT/MIT based preservative
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

### Applications

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab112618 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/1000 - 1/30000. Colorimetric: Use 1/1000-1/10000 dilution; Chemiluminescent: Use 1/1000-1/30000 dilution.
IHC-P		1/200 - 1/500.
ICC/IF		1/200 - 1/500.
ELISA		1/10000 - 1/100000. as primary antibody.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors