



## Anti- STAT1 (phospho S727) antibody [EPR3146] - (ab109461)

Application	ChIC/CUT&RUN-seq
Target	STAT1 (phospho S727)
Sample	HeLa
Treatment	IFN gamma (50ng/ml 1h)
Antibody amount	5µg
Cell number	2.5X10 <sup>5</sup>
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	hg38
Reads	10 million


### Chromosome 5



### Chromosome 1



### Chromosome 12



### Chromosome 6

