


Anti-Histone H3 (tri methyl K27) antibody [EPR18607] - ChIP Grade (ab192985)


Target	Histone H3 (tri methyl K27)
Sample	HeLa
Treatment	None
Antibody amount	4 µg
Cell number	1×10 <sup>7</sup>
Sequencing	NovaSeq 6000 (PE 150)
Peak calling	MACS 2.1.0
Reference genome	hg38

63,700 kb      63,800 kb      63,900 kb      64,000 kb      64,100 kb      64,200 kb      64,300 kb



Copyright (c) 2020 Abcam plc.

48,575 kb      48,600 kb      48,625 kb      48,650 kb



Copyright (c) 2020 Abcam plc.

