




Anti-TATA binding protein TBP antibody [EPR21954] - ChIP Grade (ab220788)

| | |
|------------------|------------------------|
| Target | TBP |
| Sample | HeLa |
| Treatment | None |
| Antibody amount | 4 µg |
| Chromatin amount | 30 µg |
| Lot no. | GR3236211 |
| Sequencing | NextSeq500, SE (75 bp) |
| Peak calling | MACS 2.1.0 |
| Reference genome | hg38 |





Chromosome 11

Input

ab220788

Refseq Genes



Chromosome 4

Input

ab220788

Refseq Genes

