



Anti-SP1 antibody [EPR22648-50] - ChIP Grade (ab231778)

Target	SP1
Sample	HeLa
Treatment	None
Antibody amount	4 µg
Chromatin amount	30 µg
Lot no.	GR3276730
Sequencing	NextSeq500, SE (75 bp)
Peak calling	MACS 2.1.0
Reference genome	hg38



7,670 kb

7,680 kb

7,690 kb


7,700 kb

Chromosome 17

Input

ab231778

Refseq Genes



229,430 kb

229,440 kb

229,450 kb

229,460 kb

229,470 kb

229,480 kb

229,490 kb

Chromosome 1

Input

ab231778

Refseq Genes

