



Anti-RNA polymerase II CTD repeat YSPTSPS (phospho S2) antibody [EPR18855-87] -  
ChIP Grade (ab238146)

Target	RNA Pol II CTD (pS2)
Sample	HeLa
Treatment	None
Antibody amount	4 µg
Chromatin amount	30 µg
Lot no.	GR3281311
Sequencing	NextSeq500, SE (75 bp)
Peak calling	MACS 2.1.0
Reference genome	hg38




kb 49,185 kb 49,190 kb 49,195 kb 49,200 kb 49,205 kb

## Chromosome 12

Input

ab238146

Refseq Genes



17,720 kb

17,730 kb

17,740 kb

## Chromosome 11

Input

ab238146

Refseq Genes



MYOD1



KCNC1

KCNC1