



## Anti- RUNX1 / AML1 antibody [EPR23044-100] - (ab240639)

|                  |                      |
|------------------|----------------------|
| Application      | ChIC/CUT&RUN-seq     |
| Target           | RUNX1 / AML1         |
| Sample           | K-562                |
| Treatment        | None                 |
| Antibody amount  | 5µg                  |
| Cell number      | 2.5X10 <sup>5</sup>  |
| Sequencing       | NovaSeq 6000 (PE150) |
| Peak calling     | MACS 2.1.0           |
| Reference genome | hg38                 |
| Reads            | 10 million           |


### Chromosome 3



### Chromosome 12




Chromosome 11



Copyright (c) 2023 Abcam plc.

Chromosome 10



Copyright (c) 2023 Abcam plc.