


Anti-Histone H3 (mutated K36M) antibody [EPR23614-91] - (ab256384)


Application	ChIC/CUT&RUN-seq
Target	Histone H3 (mutated K36M)
Sample	293T
Treatment	Transfected with H3.3K36M (Mutant) or H3.3 (WT) plasmid
Antibody amount	3µg
Cell number	2.5X10 <sup>5</sup>
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	hg38
Reads	10 million

Chromosome 3




Copyright (c) 2023 Abcam plc.

Chromosome 12




Copyright (c) 2023 Abcam plc.

Chromosome 1



Copyright (c) 2023 Abcam plc.

Chromosome 16



Copyright (c) 2023 Abcam plc.