



Anti-HNF1 alpha antibody [EPR23054-142] - (ab242140)

Application	ChIC/CUT&RUN-seq
Target	HNF1 alpha
Sample	HepG2
Treatment	None
Antibody amount	5µg
Cell number	2.5X10 ⁵
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	hg38
Reads	10 million


Chromosome 20



Chromosome 20




Chromosome 10




Copyright (c) 2023 Abcam plc.

Chromosome 16



Copyright (c) 2023 Abcam plc.

Chromosome 20



Copyright (c) 2023 Abcam plc.