


## Anti- STAT3(phosphor S727) antibody [EPR121-31] - (ab32143)


Application	ChIC/CUT&RUN-seq
Target	STAT3 (phospho S727)
Sample	HepG2
Treatment	starved overnight +IL-6(100ng/ml 30min)
Antibody amount	5µg
Cell number	2.5X10 <sup>5</sup>
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	hg38
Reads	10 million

### Chromosome 17




Copyright (c) 2023 Abcam plc.

### Chromosome 5



### Chromosome 14



### Chromosome 2

