



Anti-Histone H3 (acetyl K18) antibody [EP959Y] - ChIP Grade (ab40888)

Target	Histone H3 (acetyl K18)
Sample	HeLa
Treatment	None
Antibody amount	4 µg
Cell number	1×10 ⁷
Sequencing	NovaSeq 6000 (PE 150)
Peak calling	MACS 2.1.0
Reference genome	hg38



5,500 kb

5,525 kb

5,550 kb

5,575 kb

Input

ab40888

Refseq Genes

FBXL18

LOC221946

ACTB

FBXL18

17,690 kb

17,700 kb

17,710 kb

17,720 kb

17,730 kb

17,740 kb

17,750 kb

Input

ab40888

Refseq Genes

LINC02729

MYOD1

KCNC1

KCNC1

Copyright (c) 2021 Abcam plc.

73,400 kb

73,450 kb

73,500 kb

73,550 kb

Input

ab40888

Refseq Genes

ALB

AFP

AFP

AFM

AFM

LINC02499

R/

R/

Copyright (c) 2021 Abcam plc.

94,325 kb

94,350 kb

94,375 kb

94,400 kb

94,425 kb

Input

ab40888

Refseq Genes

SERPINA6

SERPINA2

SERPINA1

SERPINA1