



Anti-RBBP4 antibody [EPR3411] - CHIP Grade (ab79416)

Target	RBBP4
Sample	HeLa
Treatment	None
Antibody amount	3 µg
Chromatin amount	30 µg
Sequencing	NextSeq500, SE (75 bp)
Peak calling	MACS 2.1.0
Reference genome	hg38



176,075 kb

176,100 kb

176,125 kb

176,150 kb

176,175 kb

176,200 kb

Chromosome 2

Input



ab79416



Refseq Genes



Copyright (c) 2021 Abcam plc.

97,900 kb

98,000 kb

98,100 kb

98,200 kb

98,300 kb

98,400 kb

98,500 kb

Chromosome 8


Input



ab79416



Refseq Genes



Copyright (c) 2021 Abcam plc.

5,520 kb

5,530 kb

5,540 kb

5,550 kb

5,560 kb

Chromosome 7

Input



ab79416



Refseq Genes



Copyright (c) 2021 Abcam plc.