


Recombinant Rabbit IgG, monoclonal [EPR25A] - Isotype Control (ab172730)


Target	Rabbit IgG
Sample	HeLa
Treatment	None
Antibody amount	0.5ug/1ug/2ug
Cell number	10 ⁵
Sequencing	NovaSeq 6000 (PE 150)
Peak calling	MACS 2.1.0
Reference genome	hg38

Chromosome 7




Copyright (c) 2021 Abcam plc.

Chromosome 1




Copyright (c) 2021 Abcam plc.

Chromosome 6



Copyright (c) 2021 Abcam plc.

Chromosome 7



Copyright (c) 2021 Abcam plc.