


Anti-EHMT2/G9A antibody [EPR18894] (ab185050)


|                  |                       |
|------------------|-----------------------|
| Target           | EHMT2/G9A             |
| Sample           | HeLa                  |
| Treatment        | None                  |
| Antibody amount  | 8 µg                  |
| Cell number      | 1×10 <sup>7</sup>     |
| Sequencing       | NovaSeq 6000 (PE 150) |
| Peak calling     | MACS 2.1.0            |
| Reference genome | hg38                  |

Chromosome 19 (active region)




Copyright (c) 2022 Abcam plc.

Chromosome 16 (active region)




Copyright (c) 2022 Abcam plc.

### Chromosome 15 (active region)




Copyright (c) 2022 Abcam plc.

### Chromosome 4 (inactive region)



Copyright (c) 2022 Abcam plc.

### Chromosome 14 (inactive region)



Copyright (c) 2022 Abcam plc.