


Anti- Smad2+Smad3 antibody [EPR19557-4] - (ab202445)


Application	ChIC/CUT&RUN-seq
Target	Smad2+Smad3
Sample	HaCaT
Treatment	TGF-β1(7ng/ml 1h)
Antibody amount	5μg
Cell number	2.5X10 ⁵
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	hg38
Reads	10 million

Chromosome 19




Copyright (c) 2023 Abcam plc.


Chromosome 20




Chromosome 16



Chromosome 6



Chromosome 10



Copyright (c) 2023 Abcam plc.