


## Anti-ATF3 antibody [EPR22610-19] - (ab254268)


Application	ChIC/CUT&RUN-seq
Target	ATF3
Sample	Human ATF3 knockout HeLa cell line (ab264908) Human wild-type HeLa cell line (ab255448)
Treatment	None
Antibody amount	5µg
Cell number	2.5X10 <sup>5</sup>
Antibody used for Assay Quality Control	Anti-CTCF (ab188408)
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	hg38
Reads	10 million

### Chromosome 12 zoom out




Copyright (c) 2023 Abcam plc.

ab254268 Chromosome 1




Copyright (c) 2023 Abcam plc.

ab254268 Chromosome 6




Copyright (c) 2023 Abcam plc.

ab254268 Chromosome 1



Copyright (c) 2023 Abcam plc.

ab254268 Chromosome 3



Copyright (c) 2023 Abcam plc.