



Anti-Histone H3 (tri methyl K36) antibody [Y23525-232] (ab282572)

| | |
|------------------|-----------------------|
| Target | H3K36me3 |
| Sample | HeLa |
| Treatment | None |
| Antibody amount | 4 µg |
| Cell number | 1×10 ⁷ |
| Sequencing | NovaSeq 6000 (PE 150) |
| Peak calling | MACS 2.1.0 |
| Reference genome | hg38 |


Chromosome 12



Chromosome 7



Chromosome 19



Copyright (c) 2021 Abcam plc.