



Anti- Histone H3 (di methyl K4) antibody [Y47] - (ab32356)

Application	ChIC/CUT&RUN-seq
Target	Histone H3 (di methyl K4)
Sample	HeLa
Treatment	None
Antibody amount	5µg
Cell number	2.5X10 ⁵
Sequencing	NovaSeq 6000 (PE150)
Peak calling	MACS 2.1.0
Reference genome	hg38
Reads	10 million


Chromosome 16



Chromosome 16




Chromosome 7



Copyright (c) 2023 Abcam plc.

Chromosome 12



Copyright (c) 2023 Abcam plc.