

Product datasheet

Anti-PTPIP51 antibody ab69337

2 Images

Overview

Product name	Anti-PTPIP51 antibody
Description	Mouse polyclonal to PTPIP51
Host species	Mouse
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Full length protein corresponding to Human PTPIP51. Database link: NP_060615
Positive control	Human liver tissue lysate and FAM82C transfected 293T lysate.
General notes	This product was previously labelled as FAM82A2

Reproducibility is key to advancing scientific discovery and accelerating scientists' next breakthrough.

Abcam is leading the way with our range of recombinant antibodies, knockout-validated antibodies and knockout cell lines, all of which support improved reproducibility.

We are also planning to innovate the way in which we present recommended applications and species on our product datasheets, so that only applications & species that have been tested in our own labs, our suppliers or by selected trusted collaborators are covered by our Abpromise™ guarantee.

In preparation for this, we have started to update the applications & species that this product is Abpromise guaranteed for.

We are also updating the applications & species that this product has been “predicted to work with,” however this information is not covered by our Abpromise guarantee.

Applications & species from publications and Abreviews that have not been tested in our own labs or in those of our suppliers are not covered by the Abpromise guarantee.

Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, as well as customer reviews and Q&As.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	pH: 7.4 Constituent: 2.68% PBS
Purity	Protein G purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab69337** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/500 - 1/1000. Detects a band of approximately 53 kDa (predicted molecular weight: 53 kDa).

Target

Relevance	FAM82A2, or family with sequence similarity 82, member C, may participate in differentiation and apoptosis of keratinocytes. Overexpression of FAM82A2 induces apoptosis.
Cellular localization	Mitochondrion membrane; Single-pass membrane protein. Cytoplasm. Nucleus. Note=In interphase localizes in the cytoplasm, and during mitosis localizes to the spindle microtubules and spindle poles.

Images



Anti-PTPIP51 antibody (ab69337) at 1/500 dilution + Human liver tissue lysate at 25 µg


Secondary

Goat Anti-Mouse IgG (H&L)-HRP at 1/2500 dilution

Predicted band size: 53 kDa

Observed band size: 53 kDa

Additional bands at: 36 kDa. We are unsure as to the identity of these extra bands.



Western blot - Anti-PTPIP51 antibody (ab69337)

All lanes : Anti-PTPIP51 antibody (ab69337) at 1/500 dilution

Lane 1 : PTPIP51 transfected 293T lysate

Lane 2 : Non-transfected 293T lysate

Lysates/proteins at 25 µg per lane.

Secondary

Lane 1 : Non-transfected at 1/2500 dilution

Lane 2 : Non-transfected at 1/2500 dilution

Predicted band size: 53 kDa

Observed band size: 65 kDa

[why is the actual band size different from the predicted?](#)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors