

Product datasheet

Anti-PTPN22 antibody ab243320

1 Image

Overview

Product name	Anti-PTPN22 antibody
Description	Rabbit polyclonal to PTPN22
Host species	Rabbit
Tested applications	Suitable for: IHC-P
Species reactivity	Reacts with: Human
Immunogen	Recombinant fragment corresponding to Human PTPN22 aa 359-449. Sequence: SEISAKEELVLHPAKSSTSFDFLELNYSFDKNADTTMKWQ TKAFPVGEF LQKHQSLDLGSLLFEGCSNSKPVNAAGRYFNSKVPITRTK S

Database link: [Q9Y2R2](#)

 [Run BLAST with](#)

 [Run BLAST with](#)

Positive control IHC-P: Human lymph node tissue.

General notes

Reproducibility is key to advancing scientific discovery and accelerating scientists' next breakthrough.

Abcam is leading the way with our range of recombinant antibodies, knockout-validated antibodies and knockout cell lines, all of which support improved reproducibility.

We are also planning to innovate the way in which we present recommended applications and species on our product datasheets, so that only applications & species that have been tested in our own labs, our suppliers or by selected trusted collaborators are covered by our Abpromise™ guarantee.

In preparation for this, we have started to update the applications & species that this product is Abpromise guaranteed for.

We are also updating the applications & species that this product has been “predicted to work with,” however this information is not covered by our Abpromise guarantee.

Applications & species from publications and Abreviews that have not been tested in our own labs or in those of our suppliers are not covered by the Abpromise guarantee.

Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, as well as

customer reviews and Q&As.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 40% Glycerol (glycerin, glycerine), PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab243320** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		1/20 - 1/50. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

Target

Function	Seems to act on Cbl. May play a role in regulating the function of Cbl and its associated protein kinases. Acts as negative regulator of T cell receptor (TCR) signaling. Dephosphorylates and inactivates the SRC family kinases.
Tissue specificity	Predominantly expressed in lymphoid tissues and cells. Isoform 1 is expressed in thymocytes and both mature B and T-cells.
Sequence similarities	Belongs to the protein-tyrosine phosphatase family. Non-receptor class 4 subfamily. Contains 1 tyrosine-protein phosphatase domain.
Post-translational modifications	Phosphorylation on Ser-35 by PKC/PRKCD abrogates its ability to dephosphorylate and inactivate the SRC family kinases.
Cellular localization	Cytoplasm.

Images



Formalin-fixed, paraffin-embedded human lymph node tissue stained for PTPN22 using ab243320 at 1/20 dilution in immunohistochemical analysis.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-PTPN22 antibody (ab243320)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors