




Product datasheet

Anti-RAB30 antibody ab221073

1 Image

Overview

<b>Product name</b>	Anti-RAB30 antibody
<b>Description</b>	Rabbit polyclonal to RAB30
<b>Host species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> ICC/IF
<b>Species reactivity</b>	<b>Reacts with:</b> Human <b>Predicted to work with:</b> Mouse, Cow 
<b>Immunogen</b>	Recombinant fragment corresponding to Human RAB30 aa 156-202. Sequence: SDNVEKLF <del>LDLACRLISEARQNTLVNNVSSPLP</del> GEGK SISYL <del>TCCNF</del>  Database link: <a href="#">Q15771</a>   <a href="#">Run BLAST with</a>  <a href="#">Run BLAST with</a>
<b>Positive control</b>	ICC/IF - SiHa cells.

Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	pH: 7.2 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
<b>Purity</b>	Immunogen affinity purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab221073** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
ICC/IF		Use a concentration of 1 - 4 µg/ml. Recommended conditions: Fixation/Permeabilization: PFA/Triton X-100

## Target

**Sequence similarities** Belongs to the small GTPase superfamily. Rab family.

**Cellular localization** Cell membrane.

## Images



Immunofluorescent analysis of PFA fixed, Triton X-100 permeabilized SiHa cells labeling RAB30 with ab221073 at 4 µg/ml (green).

Immunocytochemistry/ Immunofluorescence - Anti-RAB30 antibody (ab221073)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.com/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors