

Product datasheet

Anti-RAB36 antibody ab191531

1 Image

Overview

Product name	Anti-RAB36 antibody
Description	Rabbit polyclonal to RAB36
Host species	Rabbit
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Mouse, Rat, Human
Immunogen	Synthetic peptide within Human RAB36. The exact sequence is proprietary. Database link: O95755
Positive control	HEK293T, NIH 3T3 and PC12 whole cell lysates.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.01% Sodium azide Constituents: 50% Glycerol, 49% PBS
Purity	Immunogen affinity purified
Purification notes	ab191531 was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen, and the purity is > 95% (by SDS-PAGE).
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab191531** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use at an assay dependent concentration. Predicted molecular weight: 36 kDa.

Target

Function	Protein transport. Probably involved in vesicular traffic.
Tissue specificity	Ubiquitously present in all tissues examined.
Sequence similarities	Belongs to the small GTPase superfamily. Rab family.
Cellular localization	Golgi apparatus membrane.

Images



All lanes : Anti-RAB36 antibody (ab191531)

Lane 1 : HEK293T whole cell lysate

Lane 2 : NIH3T3 whole cell lysate

Lane 3 : PC12 whole cell lysate

Predicted band size: 36 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors