

## Product datasheet

# Rabbit Anti-Cow IgG H&L (Texas Red ®) ab6925

### Overview

<b>Product name</b>	Rabbit Anti-Cow IgG H&L (Texas Red ®)
<b>Host species</b>	Rabbit
<b>Target species</b>	Cow
<b>Tested applications</b>	<b>Suitable for:</b> ICC/IF, IHC-P, IHC-Fr, Immunomicroscopy, Flow Cyt, ELISA
<b>Immunogen</b>	Full length native bovine IgG (purified).
<b>Conjugation</b>	Texas Red ®. Ex: 596nm, Em: 620nm

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.87% Sodium chloride, 1% BSA
<b>Purity</b>	Affinity purified
<b>Purification notes</b>	This product was prepared from monospecific antiserum by immunoaffinity chromatography using bovine IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities
<b>Conjugation notes</b>	Texas Red™ Sulfonyl Chloride (Molecular Weight 625 daltons) Absorption Wavelength: 596 nm Emission Wavelength: 620 nm Fluorochrome/Protein Ratio: 4.9 moles Texas Redä per mole of Rabbit IgG
<b>Clonality</b>	Polyclonal

### Applications

**The Abpromise guarantee** Our [Abpromise guarantee](#) covers the use of ab6925 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		

Application	Abreviews	Notes
IHC-P		
IHC-Fr		
Immunomicroscopy		
Flow Cyt		
ELISA		

#### Application notes

Flow Cyt: Use at an assay dependant concentration.

IHC-P: Use at an assay dependant concentration.

IHC-Fr: Use at an assay dependant concentration.

IM: Use at an assay dependant concentration.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors