

## Product datasheet

# Rabbit F(ab')<sub>2</sub> Anti-Cow IgG H&L (HRP) ab136757

### Overview

<b>Product name</b>	Rabbit F(ab') <sub>2</sub> Anti-Cow IgG H&L (HRP)
<b>Description</b>	Rabbit F(ab') <sub>2</sub> polyclonal Secondary Antibody to Cow IgG - H&L (HRP)
<b>Host species</b>	Rabbit
<b>Target species</b>	Cow
<b>Tested applications</b>	<b>Suitable for:</b> Flow Cyt, IHC-P, IHC-Fr, ELISA, ICC, WB
<b>Immunogen</b>	Highly purified whole cow immunoglobulin
<b>Conjugation</b>	HRP

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	Constituents: 50% Tris buffered saline, 50% Glycerol
<b>Purity</b>	Immunogen affinity purified
<b>Purification notes</b>	ab136757 is purified by affinity chromatography and is thus particularly specific.
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

### Applications

**The Abpromise guarantee** Our [Abpromise guarantee](#) covers the use of ab136757 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		Use at an assay dependent concentration.
IHC-P		1/20 - 1/100.
IHC-Fr		1/20 - 1/100.

Application	Abreviews	Notes
ELISA		1/500 - 1/5000.
ICC		1/20 - 1/100.
WB		1/100 - 1/1000.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors