

Product datasheet

Rabbit F(ab')₂ Anti-Dog IgG H&L (FITC) ab136754

Overview

Product name	Rabbit F(ab') ₂ Anti-Dog IgG H&L (FITC)
Description	Rabbit F(ab') ₂ polyclonal Secondary Antibody to Dog IgG - H&L (FITC)
Host species	Rabbit
Target species	Dog
Tested applications	Suitable for: IHC-Fr, IHC-P, Flow Cyt, ICC/IF, ICC
Immunogen	Highly purified whole immunoglobulin from dog
Conjugation	FITC. Ex: 493nm, Em: 528nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.1% Sodium azide Constituents: 99.6% PBS, 0.3% BSA
Purity	Immunogen affinity purified
Purification notes	ab136754 is purified by affinity chromatography and is thus particularly specific.
Conjugation notes	The F/P (fluorochrome/protein) molar ratio is set between 3 and 7.
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab136754 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-Fr		1/10 - 1/100.
IHC-P		1/10 - 1/100.

Application	Abreviews	Notes
Flow Cyt		1/10 - 1/100.
ICC/IF		1/50 - 1/200.
ICC		1/10 - 1/100.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors