

## Product datasheet

# Rabbit F(ab')<sub>2</sub> Anti-Sheep IgG H&L (Biotin) ab136770

### Overview

<b>Product name</b>	Rabbit F(ab') <sub>2</sub> Anti-Sheep IgG H&L (Biotin)
<b>Description</b>	Rabbit F(ab') <sub>2</sub> polyclonal Secondary Antibody to Sheep IgG - H&L (Biotin)
<b>Host species</b>	Rabbit
<b>Target species</b>	Sheep
<b>Tested applications</b>	<b>Suitable for:</b> ELISA, WB, ICC, IHC-P, IHC-Fr, Flow Cyt
<b>Immunogen</b>	Highly purified whole Sheep immunoglobulin.
<b>Conjugation</b>	Biotin

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	Constituents: 50% Glycerol, 50% PBS
<b>Purity</b>	Immunogen affinity purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

### Applications

Our [Abpromise guarantee](#) covers the use of **ab136770** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		1/1000 - 1/10000.
WB		1/200 - 1/2000.
ICC		1/50 - 1/200.
IHC-P		1/50 - 1/200.

Application	Abreviews	Notes
IHC-Fr		1/50 - 1/200.
Flow Cyt		Use at an assay dependent concentration.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors