

Product datasheet

Rabbit Anti-Goat IgG H&L (15nm Gold) ab27247

1 References

Overview

Product name	Rabbit Anti-Goat IgG H&L (15nm Gold)
Host species	Rabbit
Target species	Goat
Tested applications	Suitable for: Electron Microscopy
Immunogen	Full length Goat IgG (H&L)
Conjugation	Gold 15nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 8.20 Preservative: 0.1% Sodium azide Constituents: Tris buffered saline, 1% BSA, 30% Glycerol
Purity	Immunogen affinity purified
Clonality	Polyclonal

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of **ab27247** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Electron Microscopy		

Application notes Electron Microscopy: Use at an assay dependent dilution.

Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors