

Product datasheet

Rabbit Anti-Goat IgG+IgM+IgA H&L (HRP) ab8515

1 References

Overview

Product name	Rabbit Anti-Goat IgG+IgM+IgA H&L (HRP)
Host species	Rabbit
Target species	Goat
Specificity	Antibodies were raised in rabbit to goat IgG, IgM and IgA
Tested applications	Suitable for: ELISA, Dot blot
Conjugation	HRP

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 7.50 Constituents: 0.164% Sodium phosphate, 50% Glycerol, 0.58% Sodium chloride
Purity	Immunogen affinity purified
Purification notes	Specific antibodies were absorbed by incubation with Sepharose-bound goat immunoglobulins. Specific antibodies were eluted by acidic buffer at pH 2.5 followed by neutralisation and dialysis. After repeated binding with immobilized goat immunoglobulins a minimum of 65% protein bound.
Conjugation notes	The purified polyclonal was conjugated to horseradish peroxidase according to the periodate method followed by gel-filtration and ion-exchange chromatography to clear unbound conjugate.
Clonality	Polyclonal

Applications

Our [Abpromise guarantee](#) covers the use of **ab8515** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		1/1000.
Dot blot		1/100.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors