

Product datasheet

Rabbit Anti-Guinea pig IgG H&L ab6698

2 References

Overview

Product name	Rabbit Anti-Guinea pig IgG H&L
Host species	Rabbit
Target species	Guinea pig
Tested applications	Suitable for: IP, Immunodiffusion, Conjugation, ICC/IF, ELISA, WB, IHC-P, IHC-Fr
Immunogen	Guinea Pig IgG whole molecule

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.01% Sodium azide Constituents: 0.87% Sodium chloride, 0.164% Sodium phosphate
Purity	Affinity purified
Purification notes	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Guinea Pig IgG coupled to agarose beads.
Clonality	Polyclonal

Applications

The Abpromise guarantee Our [Abpromise guarantee](#) covers the use of ab6698 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IP		Use at an assay dependent concentration.
Immunodiffusion		Use at an assay dependent concentration.
Conjugation		Use at an assay dependent concentration.

Application	Abreviews	Notes
ICC/IF		1/1000 - 1/5000.
ELISA		1/20000 - 1/100000.
WB		1/2000 - 1/10000.
IHC-P		1/1000 - 1/5000.
IHC-Fr		1/1000 - 1/5000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors