

## Product datasheet

# Rabbit Anti-Human IgA alpha chain ab97216

★★★★★ 1 Abreviews

### Overview

<b>Product name</b>	Rabbit Anti-Human IgA alpha chain
<b>Host species</b>	Rabbit
<b>Target species</b>	Human
<b>Specificity</b>	By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgA. Cross reactivity with other immunoglobulins and light chains is less than 0.1%.
<b>Tested applications</b>	<b>Suitable for:</b> IHC-P, ELISA, WB, ICC

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	pH: 6.8 Preservative: 0.09% Sodium azide Constituent: PBS
<b>Purity</b>	Immunogen affinity purified
<b>Purification notes</b>	Antiserum was solid phase adsorbed to ensure class specificity. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads.
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

### Applications

Our [Abpromise guarantee](#) covers the use of **ab97216** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		1/200 - 1/2000.
ELISA		1/100 - 1/30000. Coating: 1/100 - 1/500; Primary: 1/1000 - 1/30000

Application	Abreviews	Notes
WB		1/1000 - 1/30000. Colorimetric: 1/1000 - 1/10000; Chemiluminescent: 1/1000 - 1/30000
ICC		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors