

Product datasheet

Rabbit Anti-Monkey IgG H&L ab112776

Overview

Product name	Rabbit Anti-Monkey IgG H&L
Description	Rabbit polyclonal Secondary Antibody to Monkey IgG - H&L
Host species	Rabbit
Target species	Monkey
Specificity	this antibody reacts specifically with monkey IgG and with light chains common to other monkey immunoglobulins. Species detected are Rhesus, Baboon, Carb Green, Cyno, Chimp and African Green. Reaction to Human was also detected. No antibody was detected against non-immunoglobulin serum proteins.
Tested applications	Suitable for: ICC/IF, WB, ELISA, IHC-P

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.20 Preservative: 0.1% Sodium azide Constituent: 99% PBS
Purity	Immunogen affinity purified
Purification notes	ab112776 was isolated by affinity chromatography using antigen coupled to agarose beads.
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab112776** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		1/200 - 1/2000.

Application	Abreviews	Notes
WB		1/1000 - 1/30000. Colorimetric: 1/1000-1/10000. Chemiluminescent: 1/1000-1/30000.
ELISA		1/100 - 1/30000. Coating: 1/100-1/500. Primary: 1/1000-1/30000.
IHC-P		1/200 - 1/2000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors