

Product datasheet

Rabbit Anti-Mouse IgG H&L (Biotin) preadsorbed ab98784

★☆☆☆☆ 1 Abreviews 1 Image

Overview

Product name	Rabbit Anti-Mouse IgG H&L (Biotin) preadsorbed
Description	Rabbit polyclonal Secondary Antibody to Mouse IgG - H&L (Biotin), pre-adsorbed
Host species	Rabbit
Target species	Mouse
Tested applications	Suitable for: IHC-P, ICC/IF, ELISA, WB
Minimal cross-reactivity	Human, Rat more details
Conjugation	Biotin

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: PBS, 0.2% BSA
Purity	Immunogen affinity purified
Purification notes	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Biotin.
Clonality	Polyclonal
Isotype	IgG


Applications

Our [Abpromise guarantee](#) covers the use of **ab98784** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		1/500 - 1/5000.
ICC/IF		1/200 - 1/1000.
ELISA		1/20000 - 1/400000. as primary
WB		1/20000 - 1/400000. chemiluminescent detection

Images



Immunohistochemical analysis of dog skin tissue, detecting staining of NQO1 with ab98784.

Tissue was fixed with formaldehyde and blocked with none for ; antigen retrieval was by heat mediation in a citrate buffer (pH 6). Samples were incubated with primary antibody ([ab28947 Mouse monoclonal \[A180\] to NQO1](#) at 1/175 in BSA in TBS) for 45 minutes. ab98784 (1/500) was used as the secondary antibody.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Rabbit Anti-Mouse IgG H&L (Biotin) preadsorbed (ab98784)

This image is courtesy of an anonymous Abreview

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors