

Product datasheet

RBM20 peptide ab126284

Description

Product name	RBM20 peptide
Animal free	No
Nature	Synthetic

Specifications

Our [Abpromise guarantee](#) covers the use of **ab126284** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Blocking - Blocking peptide for Recombinant Human ITLN1 protein (denatured) (ab109850)
Form	Lyophilised

Preparation and Storage

Stability and Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Reconstitution	Resuspend in water or in neutral buffer

General Info

Tissue specificity	Expressed in the heart.
Involvement in disease	Defects in RBM20 are the cause of cardiomyopathy dilated type 1DD (CMD1DD) [MIM:613172]. A disorder characterized by ventricular dilation and impaired systolic function, resulting in congestive heart failure and arrhythmia. Patients are at risk of premature death.
Sequence similarities	Contains 1 RRM (RNA recognition motif) domain.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors