

Product datasheet

Recombinant 46 kDa surface antigen protein (His tag)
 ab240866

1 Image

Description

Product name	Recombinant 46 kDa surface antigen protein (His tag)	
Purity	> 90 % SDS-PAGE.	
Expression system	Yeast	
Accession	P0C0J7	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Sequence	CGQTESGSTSDSKPQAETLKHKVSNDSIRIALTDPDNPRW ISAQKDIISY VDETEAATSTITKNQDAQNNWLTQQANLSPAPKGFIIAPEN GSGVGTAVN TIADKGIPVAYDRLITGSDKYDWVVSFDNEKVGELQGLSLA AGLLGKED GAFDSIDQMNEYLKSHMPQETISFYTIAGSQDDNNSQYFYN GAMKVLKEL MKNSQNKIIDLSPEGENAVYVPGWNYGTAGQRIQSFLTINK DPAGGNKIK AVGSKPASIFKGFLAPNDGMAEQAITKLKLEGFDTQKIFVT GQDYNDKAK TFIKGDQNMTYKPKDKVLGKVAVEVLRVLIAKKNKASRSE VENELKAKL PNISFKYDNQTYKVQGKNINTILVSPVIVTKANVDNPDA	
Predicted molecular weight	45 kDa including tags	
Amino acids	28 to 416	
Tags	His tag N-Terminus	
Additional sequence information	Full length mature chain without signal peptide. Mycoplasma hypopneumoniae (strain 232).	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab240866** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid


Preparation and Storage

Stability and Storage Shipped at 4°C. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.2

Constituents: Tris buffer, 50% Glycerol (glycerin, glycerine)

Images



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel analysis of ab240866.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors