

Product datasheet

Recombinant Chicken TL1A protein ab90800

Description

Product name	Recombinant Chicken TL1A protein	
Purity	> 95 % SDS-PAGE. Purified by Ion-exchange chromatography	
Expression system	Yeast	
Protein length	Protein fragment	
Animal free	No	
Nature	Recombinant	
Species	Chicken	
Sequence	ERSSHFLKQR AVAAVTDTLP SAEKPRAHLT VKKQEPSSTT GSHLPILQWE DKRGLAFTKN NLSYSSNALV IPVSGDYVYVY AQVTFRGPSD TSSKTSSVTA VITKVTDSYP EPTQLLTSTK TLSEERNNWF QPMLGAVVS LEIGDKLMVN VSDIKLVDYT KEHKTFFGAF LL	
Amino acids	68 to 239	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab90800** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Lyophilized

Preparation and Storage

Stability and Storage Shipped at 4°C. Upon delivery aliquot. Upon reconstitution add a carrier protein (0.1% BSA). Store at -20°C. Avoid freeze / thaw cycle. For long term storage it is recommended to add a carrier protein on reconstitution (0.1% HSA or BSA). Stable for 12 months at -20°C. Working aliquots stored with a carrier protein are stable for at least 3 months at -20°C to -80°C..

Constituents: 10% Trehalose, 2.68% PBS

Reconstitution Reconstitute with sterile phosphate-buffered saline containing at least 0.1% carrier protein.

General Info

Function	Receptor for TNFRSF25 and TNFRSF6B. Mediates activation of NF-kappa-B. Inhibits vascular endothelial growth and angiogenesis (in vitro). Promotes activation of caspases and apoptosis.
Tissue specificity	Specifically expressed in endothelial cells. Detected in monocytes, placenta, lung, liver, kidney, skeletal muscle, pancreas, spleen, prostate, small intestine and colon.
Sequence similarities	Belongs to the tumor necrosis factor family.
Cellular localization	Secreted and Membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors