

Product datasheet

Recombinant Cytomegalovirus pp150 protein ab43039

Description

Product name	Recombinant Cytomegalovirus pp150 protein
Biological activity	Immunoreactive with sera of CMV infected individuals.
Purity	> 95 % SDS-PAGE. Sephacrose derived purification.
Expression system	Escherichia coli
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Additional sequence information	The protein contains the CMV Pp150 (UL32) immunodominant regions.

Specifications

Our [Abpromise guarantee](#) covers the use of **ab43039** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA Western blot
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.395% Tris HCl, 0.0292% EDTA, 50% Glycerol
------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------

General Info

Relevance	Cytomegalovirus (CMV), is a genus of Herpes viruses; in humans the species is known as Human herpesvirus 5 (HHV-5). It belongs to the Betaherpesvirinae subfamily of Herpesviridae. The name means "very big cell virus".
------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors