

Product datasheet

Recombinant Fe/S biogenesis protein NfuA ab225963

Description

Product name	Recombinant Fe/S biogenesis protein NfuA	
Purity	> 90 % SDS-PAGE.	
Expression system	Escherichia coli	
Accession	Q8DDU2	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Vibrio vulnificus	
Sequence	MSNITTEAAQTHFANLLGQQPDGNTNIRVFVVNPGTQNA ECGVSYCPPEA VEATDTEIPYQSFSAYVDELSLPFLEDAEIDYVTDKMG SQLTLKAPNAKM RKVADDAPLLERVEYAIQTQVNPQLAGHGGHVKLMEIT DAGVAMAFGGG CNGCSMVDVTLKEGIEKELLQQFSGELTAVRDATEHD RGDHSYY	
Predicted molecular weight	21 kDa	
Amino acids	1 to 194	
Additional sequence information	Strain CMCP6.	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab225963** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
 Constituents: 50% Glycerol, Tris buffer

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors