

## Product datasheet

# Recombinant HIV1 integrase protein ab49075

### Description

---

<b>Product name</b>	Recombinant HIV1 integrase protein
<b>Biological activity</b>	This protein reacts strongly with human HIV positive serum.
<b>Purity</b>	> 95 % SDS-PAGE.
<b>Expression system</b>	Escherichia coli
<b>Protein length</b>	Protein fragment
<b>Animal free</b>	No
<b>Nature</b>	Recombinant

### Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab49075** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	ELISA Western blot
<b>Form</b>	Liquid
<b>Additional notes</b>	This protein reacts strongly with human HIV positive serum.  Abcam has not and does not intend to apply for the REACH Authorisation of customers' uses of products that contain European Authorisation list (Annex XIV) substances. It is the responsibility of our customers to check the necessity of application of REACH Authorisation, and any other relevant authorisations, for their intended uses.

### Preparation and Storage

---

<b>Stability and Storage</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.  pH: 7.2 Constituents: 0.2% Triton-X-100, 9% Urea, 0.395% Tris HCl, 50% Glycerol (glycerin, glycerine)
------------------------------	--

### General Info

---

<b>Relevance</b>	Integrase is an enzyme, a product of the HIV1 Gag/Pol gene; the other two enzymes being reverse
------------------	---

transcriptase and protease. The human immunodeficiency virus (HIV-1) uses an enzyme, a so-called integrase, to carry out the integration of its viral DNA into the host chromosome thereby tricking the host cell machinery into making viral proteins. The HIV-1 integrase is a protein of 32 kDa and is composed of three domains. No cellular homologue of HIV integrase has been described, so potential inhibitors to this enzyme could be relatively nontoxic.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors