

## Product datasheet

# Recombinant HIV2 gp32 protein ab49083

### Description

---

<b>Product name</b>	Recombinant HIV2 gp32 protein
<b>Purity</b>	> 95 % SDS-PAGE.
<b>Expression system</b>	Escherichia coli
<b>Protein length</b>	Protein fragment
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	HIV-2

### Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab49083** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	ELISA Western blot
<b>Form</b>	Liquid

### Preparation and Storage

---

<b>Stability and Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 7.2 Preservative: 0.01% Sodium azide Constituents: 9% Urea, 0.395% Tris HCl, 0.0292% EDTA, 50% Glycerol (glycerin, glycerine)
------------------------------	--

### General Info

---

<b>Relevance</b>	Two species of HIV infect humans: HIV-1 and HIV-2. HIV-1 is more virulent. It is easily transmitted and is the cause of the majority of HIV infections globally. HIV-2 is less transmittable and is largely confined to West Africa. HIV1 and HIV2 have similar gag (viral core) and pol (polymerase) regions, they have relatively dissimilar env (envelope) regions. Owing to this lack of homology in the envelope region, there is little serologic cross-reactivity of the antibodies directed against the envelope antigens of both HIV1 and HIV2. gp 32 is a form of truncated gp41; gp41 is formed from
------------------	---

the envelope protein, gp160.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors