

Product datasheet

Recombinant HPV16 E6 protein (His tag) ab226447

1 References 3 Images

Description

| | | |
|--|---|--|
| Product name | Recombinant HPV16 E6 protein (His tag) | |
| Purity | > 90 % SDS-PAGE. | |
| Expression system | Escherichia coli | |
| Accession | P03126 | |
| Protein length | Full length protein | |
| Animal free | No | |
| Nature | Recombinant | |
| Sequence | MHQKRTAMFQDPQERPRKLPQLCTELQTTIHDILECVYCK QLLRREVY DFAFRDLCIVYRDGNPYAVCDKCLKFYKISEYRHYCYSLY GTTLEQQYN KPLCDLLIRCINCQKPLCPEEKQRHLDKKQRFHNIRGRWT GRCMCCRSS RTRRETQL | |
| Predicted molecular weight | 23 kDa including tags | |
| Amino acids | 1 to 158 | |
| Tags | His tag N-Terminus | |
| Additional sequence information | Human papillomavirus type 16 | |

Specifications

Our [Abpromise guarantee](#) covers the use of **ab226447** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|--------------------------|-------------------|
| Applications | SDS-PAGE |
| | Mass Spectrometry |
| Mass spectrometry | LC-MS/MS |
| Form | Liquid |

Preparation and Storage

Stability and Storage Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors