

Product datasheet

Recombinant Human ALKBH7 protein ab164545

1 Image

Description

<b>Product name</b>	Recombinant Human ALKBH7 protein	
<b>Expression system</b>	Wheat germ	
<b>Protein length</b>	Full length protein	
<b>Animal free</b>	No	
<b>Nature</b>	Recombinant	
<b>Species</b>	Human	
<b>Sequence</b>	MAGTGLLALRTLPGPSWVRGSGPSVLSRLQDAAVVRPG FLSTAEETLSR ELEPELRRRRYEYDHWDAAIHGFRETEKSRWSEASRAILQ RVQAAAFGPG QTLSSVHVLDLEARGYIKPHVDSIKFCGATIAGLSLLSPSV MRLVHTQE PGEWLELLLEPGSLYLGRSARYDFSHEILRDEESFFGERR IPRGRRISV ICRSLPEGMGPGESGQPPAC	
<b>Amino acids</b>	1 to 221	
<b>Tags</b>	GST tag N-Terminus	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab164545** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	Western blot
	ELISA
<b>Form</b>	Liquid

Additional notes

Preparation and Storage

<b>Stability and Storage</b>	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00
------------------------------	---

Constituents: 0.31% Glutathione, 0.79% Tris HCl

## General Info


---

<b>Function</b>	Probable dioxygenase that requires molecular oxygen, alpha-ketoglutarate and iron.
<b>Tissue specificity</b>	Widely expressed, with highest expression in pancreas, followed by spleen, prostate, ovary and placenta.
<b>Sequence similarities</b>	Belongs to the alkB family.
<b>Cellular localization</b>	Cytoplasm. Nucleus. Secreted. Has a predicted N-terminal signal sequence, indicating it may be secreted. Detected in cytoplasm and nucleus when expressed as fusion protein with an N-terminal tag.

---

## Images

---



ab164545 on a 12.5% SDS-PAGE stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors