

Product datasheet

Recombinant Human ANGPTL7 protein ab160886

1 Image

Overview

Product name	Recombinant Human ANGPTL7 protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Wheat germ
Amino Acid Sequence	
Species	Human
Sequence	MLKKPLSAVTWLCIFIVAFVSHPAWLQKLSKHKTPAQP QLKAANCCEEVK ELKAQVANLSSLLSELNKKQERDWVSVMQVMELES NSKRME SRLTDAES KYSEMNNQIDIMQLQAAQTVTQTSADAIYDCSSLYQKN YRISGVYKLPPD DFLGSPELEVFCDMETSGGGWTIIQRRKSGLVSFYRD WKQYKQGFSGIRG DFWLGNEHIHRLSRQPTRLRVEMEDWEGNLRYAEYSH FVLGNELNSYRLF LGNYTGNVGN DALQYHNNTAFSTKDKDNDNCLDKCA QLRKGGYWYNCCTD SNLNGVYYRLGEHNKHLDGITWYGWHGSTYSLKRVEM KIRPEDFKP
Amino acids	1 to 346
Tags	GST tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab160886** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA
	Western blot

Form	Liquid
-------------	--------

Additional notes

Protein concentration is above or equal to 0.05 mg/ml.

Preparation and Storage

Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Tissue specificity

Highly and specifically expressed in the cornea where is confined to the stromal layer.


Sequence similarities

Contains 1 fibrinogen C-terminal domain.

Cellular localization

Secreted.

Images



ab160886 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human ANGPTL7
protein (ab160886)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors