

Product datasheet

Recombinant Human Bag3 protein ab95384

1 Image

Description

Product name	Recombinant Human Bag3 protein
Purity	> 90 % SDS-PAGE. ab95384 is purified using conventional chromatography techniques.
Expression system	Escherichia coli
Accession	<u>O95817</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	<pre> MGSSHHHHHH SSGLVPRGSH MSAATHSPMM QVASGNGDRD PLPPGWEIKI DPQTGWPFV DHNSRTTWN DPRVPSEGPK ETPSSANGPS REGSRLPPAR EGHPVYPLR PGYIPIVLH EGAENRQVHP FHVYPQGMQ RFRTEAAAAA PQRSQSPLRG MPETTQPKQ CGQVAAAAA QPPASHGPER SQSPAASDCS SSSSSASLPS SGRSSLGSHQ LPRGYISIPV IHEQNVTRPA AQPSFHQAQK THYPAQQGEY QTHQPVYHKI QGDDWEPRPL RAASPFRSSV QGASSREGSP ARSSTPLHSP SPIRVHTVVD RPQQPMTHRE TAPVSQPENK PESKPGVGP ELPPGHIPIQ VIRKEVDSKP VSQKPPPSE KVEVKVPPAP VPCPPSPGP SAVPSSPKSV ATEERAAPST APAEATPPKP GEAEAPPKHP GVLKVEAILE KVQGLEQAVD NFEGKKTDDK YLMIEEYLT ELLALDSVDP EGRADVRQAR RDGVRKVQTI LEKLEQKAID VPGQVQVYEL QPSNLEADQP LQAIMEMGAV AADKGKKNAG NAEDPHTETQ QPEATAAATS NPSSMTDTPG NPAAP </pre>
Amino acids	1 to 575
Tags	His tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab95384** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.00174% PMSF, 0.316% Tris HCl, 0.0292% EDTA, 10% Glycerol


General Info

Function Inhibits the chaperone activity of HSP70/HSC70 by promoting substrate release. Has anti-apoptotic activity.

Involvement in disease Defects in BAG3 are the cause of myopathy myofibrillar BAG3-related (MFM-BAG3) [MIM:612954]. A neuromuscular disorder that results in early-onset, severe, progressive, diffuse muscle weakness associated with cardiomyopathy, severe respiratory insufficiency during adolescence, and a rigid spine in some patients. At ultrastructural level, muscle fibers display structural alterations consisting of replacement of the normal myofibrillar markings by small, dense granules, or larger hyaline masses, or amorphous material.

Sequence similarities Contains 1 BAG domain.
Contains 2 WW domains.

Images



15% SDS-PAGE analysis of 3µg ab95384.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors