

## Product datasheet

# Recombinant Human BCL7C protein ab174409

1 Image

### Overview

<b>Product name</b>	Recombinant Human BCL7C protein
<b>Protein length</b>	Full length protein

### Description

<b>Nature</b>	Recombinant
<b>Source</b>	Escherichia coli
<b>Amino Acid Sequence</b>	
<b>Accession</b>	<a href="#">Q8WUZO</a>
<b>Species</b>	Human
<b>Sequence</b>	<p>MGSSHHHHHHSSGLVPRGSHMGSMAGRTVRAETRS            RAKDDIKKVMATIEK            VRRWEKRWVTVGDTSLRIFKWVPVVDPQEEERRRAG            GGAERSRGRERRGR            GASPRGGGPLILLDLNDENSNQSFHSEGLQKGTEPS            PGGTPQPSRPVSP            AGPPEGVPEEAQPPRLGQERDPGGITAGSTDEPPML            TKEEPVPELLEAEA            PEAYPVFEPVPPVPEAAQGDTEDESEGAPPLKRICPNA            PDP</p>
<b>Molecular weight</b>	26 kDa including tags
<b>Amino acids</b>	1 to 217
<b>Tags</b>	His tag N-Terminus

### Specifications

Our [Abpromise guarantee](#) covers the use of **ab174409** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	<p>SDS-PAGE</p> <p>Mass Spectrometry</p>
<b>Mass spectrometry</b>	MALDI-TOF
<b>Purity</b>	> 85 % SDS-PAGE.

**Form** Liquid

## Preparation and Storage

**Stability and Storage** Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.


## General Info

**Function** May play an anti-apoptotic role.

**Tissue specificity** Ubiquitous.

**Sequence similarities** Belongs to the BCL7 family.

## Images



15% SDS-PAGE analysis of ab174409 (3µg).

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors