

Product datasheet

Recombinant Human Bradykinin protein ab182820

Description

Product name	Recombinant Human Bradykinin protein	
Purity	<p>> 90 % SDS-PAGE.</p> <p>Expressed in E.coli as inclusion bodies. Final product was refolded using a unique “temperature shift inclusion body refolding” technology and chromatographically purified.</p>	
Expression system	Escherichia coli	
Accession	P01042-2	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	<pre> MASMTGGQQMGRGHHHHHHGNLYFQGGEFQESQSEEID CNDKDLFKAVDA ALKKYNSQNSNNQFVLYRITEATKTVGSDFYFKYEIKE GDCPVQSGK TWQDCEYKDAAKAATGECTATVGKRSSTKFSVATQTCQI TPAEGPVVTAQ YDCLGCVHPISTQSPDLEPILRHGIQYFNNNTQHSSLFMLN EVKRAQRQV VAGLNFRMTYSIVQTNCSKENFLFLTPDCKSLWNGDTGE CTDNAYIDIQL RIASFQNCDIYPGKDFVQPPTKICVGCPRDIPTNSPELEE TLHTITKL NAENNATFYFKIDNVKKARVQVVAGKKYFIDFVARETTCS KESNEELTES CETKKLGGSLDCNAEVYVVPWEKKIYPTVNCQPLGMISLM K </pre>	
Predicted molecular weight	48 kDa including tags	
Amino acids	19 to 380	
Tags	His tag N-Terminus , T7 tag N-Terminus	
Additional sequence information	Corresponding to mature Isoform LMW. Constructed with codon optimization and expressed with a small T7-His-TEV cleavage site Tag (29aa) fusion at its N-terminal.	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab182820** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituent: 0.32% Tris HCl

Contains NaCl, EDTA, KCl, Arginine, DTT and Glycerol.

General Info

Function (1) Kininogens are inhibitors of thiol proteases; (2) HMW-kininogen plays an important role in blood coagulation by helping to position optimally prekallikrein and factor XI next to factor XII; (3) HMW-kininogen inhibits the thrombin- and plasmin-induced aggregation of thrombocytes; (4) the active peptide bradykinin that is released from HMW-kininogen shows a variety of physiological effects: (4A) influence in smooth muscle contraction, (4B) induction of hypotension, (4C) natriuresis and diuresis, (4D) decrease in blood glucose level, (4E) it is a mediator of inflammation and causes (4E1) increase in vascular permeability, (4E2) stimulation of nociceptors (4E3) release of other mediators of inflammation (e.g. prostaglandins), (4F) it has a cardioprotective effect (directly via bradykinin action, indirectly via endothelium-derived relaxing factor action); (5) LMW-kininogen inhibits the aggregation of thrombocytes; (6) LMW-kininogen is in contrast to HMW-kininogen not involved in blood clotting.

Tissue specificity Secreted in plasma. T-kinin is detected in malignant ovarian, colon and breast carcinomas, but not in benign tumors.

Involvement in disease Defects in KNG1 are the cause of high molecular weight kininogen deficiency (HMKW deficiency) [MIM:228960]. HMKW deficiency is an autosomal recessive coagulation defect. Patients with HMKW deficiency do not have a hemorrhagic tendency, but they exhibit abnormal surface-mediated activation of fibrinolysis.

Sequence similarities Contains 3 cystatin domains.

Post-translational modifications Bradykinin is released from kininogen by plasma kallikrein. Hydroxylation of Pro-383 occurs prior to the release of bradykinin. Phosphorylation sites are present in the extracellular medium. N- and O-glycosylated. O-glycosylated with core 1 or possibly core 8 glycans.

Cellular localization Secreted > extracellular space.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors