

Product datasheet

Recombinant Human C2orf60 protein ab139236

1 Image

Description

Product name	Recombinant Human C2orf60 protein	
Purity	> 95 % SDS-PAGE. ab139236 is purified using conventional chromatography techniques.	
Expression system	Escherichia coli	
Accession	A2RUC4	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	MGSSHHHHHHSSGLVPRGSHMGSMAGQHLPVPRLEGVS REQFMQHLYPQR KPLVLEGIDLGPCTSKWTVDYLSQVGGKKEVKIHVAAVAQ MDFISKNFVY R TL PFDQLVQRAAEEKHKEFFVSEDEKYYLRSLGEDPRK DVADIRKQFPL LKGDIKFPEFFKEEQFFSSVFRISSPGLQLWTHYDVMNDL LIQVTGKKRV VLFSPRDAQLYLKGTKSEVLNIDNPD LAKYPLFSKARRYE CSLEAGDVL FIPALWFHNVISEEFGVGVNIFWKHLPSECYDKTDTYGNK DPTAASRAAQ ILDRALKTLAELPEEYRDFYARRMVLHIQDKAYSKNSE	
Predicted molecular weight	39 kDa including tags	
Amino acids	1 to 315	
Tags	His tag N-Terminus	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab139236** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
	Mass Spectrometry

Mass spectrometry	MALDI-TOF
Form	Liquid


Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 0.02% DTT, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.88% Sodium chloride
------------------------------	--

General Info

Relevance	The function of the C2orf60 protein remains unknown.
------------------	--

Images



15% SDS-PAGE analysis of 3 µg ab139236.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors