

Product datasheet

Recombinant Human CCR8 protein ab158145

1 Image

Description

Product name	Recombinant Human CCR8 protein
Expression system	Wheat germ
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human

Sequence	<p>MDYTLDSLVTTVTDYYPDIFSSPCDAELIQTNGKLLLA VFYCLLFVFSL LGNSLVILVLVVCKKLRSDVYLLNLALSDLLFVFSFPF QTYLLDQWV FGTVMCKVVVSGFYIGFYSSMFFITLMSVDRYLAVVHA VYALKVRTIRMG TTLCLAVWLTAIMATIPLLVFYQVASEDGVLQCYSFYNQ QTLKWKIFTNF KMNILGLLIPFTIFMFCYKILHQLKRCQNHNKTKAIRLVLI VVIASLLF WVPFNVVLFSLTSLHSMHILDGCSISQQLTYATHVTEIISF THCCVNPVIY AFVGEKFKKHLSEIFQKSCSQIFNYLGRQMPRESCEK SSSCQQHSSRSSS VDYIL</p>
-----------------	---

Amino acids	1 to 355
Tags	GST tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab158145** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA
	Western blot

Form	Liquid
-------------	--------

Additional notes

Preparation and Storage

Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function

Receptor for the chemokine CCL1/SCYA1/I-309. May regulate monocyte chemotaxis and thymic cell line apoptosis. Alternative coreceptor with CD4 for HIV-1 infection.


Sequence similarities

Belongs to the G-protein coupled receptor 1 family.

Cellular localization

Cell membrane.

Images



ab158145 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human CCR8 protein
(ab158145)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors