

Product datasheet

Recombinant Human CD134 / OX40L receptor protein (Tagged) ab271642

1 Image

Description

Product name	Recombinant Human CD134 / OX40L receptor protein (Tagged)	
Purity	>= 90 % SDS-PAGE.	
Expression system	HEK 293 cells	
Accession	P43489	
Protein length	Protein fragment	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	LHCVGDYTPSNDRCCHECRPGNGMVSRCSRSQNTVCRP CGPGFYNDVVSS KPCKPCTWCNLRSGSERKQLCTATQDTVCRCRAGTQPL DSYKPGVDCAPC PPGHFSPGDNQACKPWTNCTLAGKHTLQPASNSSDAICE DRDPPATQPQE TQGPPARPITVQPTEAWPRTSQQPSTRPVEVPGGRAVA	
Predicted molecular weight	70 kDa	
Molecular weight information	This protein runs at a higher MW by SDS-PAGE due to glycosylation.	
Amino acids	29 to 216	
Tags	Fc tag C-Terminus	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab271642** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid

Preparation and Storage

Stability and Storage

Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle.

pH: 8

Constituents: 0.63% Tris HCl, 0.64% Sodium chloride, 0.02% Potassium chloride, 20% Glycerol (glycerin, glycerine)

General Info

Function

Receptor for TNFSF4/OX40L/GP34.


Sequence similarities

Contains 4 TNFR-Cys repeats.

Cellular localization

Membrane.

Images



SDS-PAGE analysis of ab271642.

This protein runs at a higher MW by SDS-PAGE due to glycosylation.

SDS-PAGE - Recombinant Human CD134 / OX40L receptor protein (Tagged) (ab271642)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors