

Product datasheet

Recombinant Human CDH11 protein ab192296

Description

Product name	Recombinant Human CDH11 protein	
Purity	> 95 % SDS-PAGE. Purity greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.	
Endotoxin level	< 1.000 Eu/μg	
Expression system	HEK 293 cells	
Accession	P55287	
Protein length	Protein fragment	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	<p>FAPERRGHLRPSFHGHHEKGKEGQVLQRSKRGWVWNQ FFVIEEYTGPDV LVGRLHSDIDSGDGNIKYILSGEGAGTIFIDDKSGNIHATKT LDREERA QYTLMAQAVDRDTNRPLEPPSEFMKVQDINDNPPEFLHE TYHANVPERS NVGTSVIQVTASDADDPTYGNSAKLVYSILEGQPYFSVEA QTGIIRTALP NMDREAKEEYHVVIQAKDMGGHMGGLSGTTKVMITLTDV NDNPPKFPQSV YQMSVSEAAVPGEEVGRVKAKDPDIGENGLVTYNIVDGD GMESFEITTDY ETQEGVIKLLKPVDFETKRAYSLKVEAANVHIDPKFISNGP FKDTVTVKI AVEDADEPPMFLAPSYIHEVQENAAAGTVVGRVHAKDPD AANSPIRYSID RHTDLDRFF TINPEDGFIKTTKPLDREETAWLNITVFAAEIH NRHQEAKV PVAIRVLDVNDNAPKFAAPYEGFICESDQTKPLSNQPVTIS ADDKDDTA NGPRFIFSLPPEIIHNPNTVRDNRDNTAGVYARRGGFSRQ KQDLYLLPI VISDGGIPPMSSTNTLTIKVCGCDVNGALLSCNAEAYILNA GLSTVDDIE GRMDEPKSCDKTHTCPPCPAPPELLGGPSVFLFPPKPKD</p>	

TLMISRTPEVTC
VVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNST
YRVVSVLTVLHQ
DWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTL
PPSREEMTKN
QVSLTCLVKGFPYSDIAVEWESNGQPENNYKTTPPVLDS
DGSFFLYSKLT
VDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSPGKHH
HHHH

Predicted molecular weight	94 kDa including tags
Amino acids	23 to 617
Tags	His tag C-Terminus
Additional sequence information	Fused with a FC-6His tag at the C-terminus.

Specifications

Our [Abpromise guarantee](#) covers the use of **ab192296** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE HPLC
Form	Lyophilized

Preparation and Storage

Stability and Storage	Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 7.40 Constituent: 100% PBS Lyophilized from a 0.2 µM filtered solution.
------------------------------	---

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors