

## Product datasheet

# Recombinant Human CLEC14A protein ab182814

### 1 References

#### Description

---

<b>Product name</b>	Recombinant Human CLEC14A protein
<b>Purity</b>	> 90 % SDS-PAGE. The final product was refolded using unique “temperature shift inclusion body refolding” technology and chromatographically purified.
<b>Expression system</b>	Escherichia coli
<b>Accession</b>	<a href="#">Q86T13</a>
<b>Protein length</b>	Protein fragment
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	Human
<b>Sequence</b>	<pre> MASMTGGQQMGRGHHHHHHENLYFQGGEFEHPTADRAG CSASGACYSLHH ATMKRQAAEEACILRGGALSTVRAGAELRAVLALLRAGPG PGGGSKDLLF WVALERRRSHCTLENEPLRGFSWLSSDPGGLESRTLQW VEEPQRSC TARR CAVLQATGGVEPAGWKEMRCHLRANGYLCKYQFEVLCP APRPGAASNLSY RAPFQLHSAALDFSPPGTEVSALCRGQLPISVTCIADEIGA RWDKLSGDV LCPCPGRYLRAGKCAELPNCLDDLGGFACECATGFELGK DGRSCVTSGEG QPTLGGTGVPTRRPPATATSPVPQRTWPIRVDEKLGETPL VPEQDNSVTS IPEIPRWGSQSTMSTLQMSLQAESKATITPSGVSISKFNST TSSATPQAF DSSSAV </pre>
<b>Predicted molecular weight</b>	43 kDa including tags
<b>Amino acids</b>	22 to 398
<b>Tags</b>	His-T7 tag N-Terminus
<b>Additional sequence information</b>	T7-His-TEV cleavage site Tag (29aa) fusion at its N-terminal. Extracellular domain. NP_778230.

#### Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab182814** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** SDS-PAGE

**Form** Liquid

---

## Preparation and Storage

**Stability and Storage** Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituent: 0.32% Tris HCl

Contains NaCl, KCl, EDTA, arginine, DTT and Glycerol.

---

## General Info

**Sequence similarities** Contains 1 C-type lectin domain.

Contains 1 EGF-like domain.

**Cellular localization** Membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

---

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

---

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors